


# GEM FLOW COUPLER

## Device and System

<b>EN</b>	Instructions for use .....	7
<b>FR</b>	Mode d'emploi.....	14
<b>DE</b>	Gebrauchsanweisung .....	21
<b>IT</b>	Istruzioni per l'uso.....	28
<b>ES</b>	Instrucciones de uso.....	35
<b>NL</b>	Gebruiksaanwijzing .....	42
<b>DA</b>	Betjeningsvejledning.....	49
<b>SV</b>	Bruksanvisning .....	56
<b>TR</b>	Kullanım Talimatları.....	63
<b>NO</b>	Bruksanvisning .....	70
<b>PT</b>	Instruções de utilização.....	77
<b>EL</b>	Οδηγίες χρήσης.....	84
<b>PL</b>	Instrukcja użycia.....	91

**A**  
FLOW COUPLER Device

**A**

FLOW COUPLER Device  
Dispositif FLOW COUPLER  
FLOW COUPLER-Gerät  
Dispositivo FLOW COUPLER  
Dispositivo FLOW COUPLER  
FLOW COUPLER-hulpmiddel  
FLOW COUPLER-enhed  
FLOW COUPLER-enhet  
FLOW COUPLER Cihazı  
FLOW COUPLER-enhet  
Dispositivo FLOW COUPLER  
Συσκευή FLOW COUPLER  
Urządzenie FLOW COUPLER

**B**

External Lead  
Fil externe  
Externes Kabel  
Elettrocattetero esterno  
Conductor eléctrico externo  
Externe kabel  
Ekstern ledning  
Extern ledning  
Harici Kablo  
Ekstern leder  
Condutor Externo  
Εξωτερική απαγωγή  
Odprowadzenie zewnętrzne

**C**

Jaw Assembly  
Assemblage de la mâchoire  
Haltevorrichtung  
Gruppo ganasce  
Ensamblaje de la mordaza  
Klemsamenstel  
Kæbemontage  
Kåftverktyg  
Ağız Takımı  
Kjevemontering

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
Figure 20

## SYMBOL GLOSSARY PER US FD&C ACT:

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
----------	--------	--------------	----------------	---------------

The following symbols and definitions pertain to the GEM FLOW COUPLER Device:

ISO-15223-1*		Consult instructions for use	Consult instructions for use	5.4.3
ISO-15223-1		Do not reuse	Do not reuse	5.4.2
ISO-15223-1		Use by Date	Use by Date	5.1.4
ISO-15223-1		Sterilized using ethylene oxide	Sterilized using ethylene oxide	5.2.3
ISO-15223-1		Do not use if package is damaged	Do not use if the product sterilization barrier or its packaging is compromised	5.2.8

The following symbols and definitions pertain to the GEM FLOW COUPLER Device and System:

ISO-15223-1		Caution	Attention, consult accompanying documents.	5.1.5
			This product and package is not made with natural rubber latex.	
			Content	
			CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.	
IEC 60417			Type CF Applied Part	5355
			RF Transmitter	
IEC 60417			Direct current	5031
ISO-15223-1		Catalog number	Catalog number	5.1.6
ISO-15223-1		Batch code	Lot number	5.1.5
ISO-15223-1		Manufacturer	Manufacturer	5.1.1
ISO-15223-1		Authorized Representative in the European Community	Authorized Representative in the European Community	5.1.2

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO 7010			Refer to instruction booklet (symbol white on blue)	M002
ISO-15223-1		Do not re-sterilize	Do not re-sterilize	5.2.6
ISO-15223-1		Date of manufacture	Date of manufacture	5.1.3

\*ISO 15223-1:2016 Medical Devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

### Additional symbols and graphics on the product labeling that are not required by the US FD&C Act:

Symbol	Symbol Description
	Made in the U.S.A.
	Manufacturer part number
	Internal tracking number

Graphic	Graphic Meaning
	Size of the GEM FLOW COUPLER Device (inner diameter of the FLOW COUPLER rings)

### DESCRIPTION:

The Synovis MCA GEM FLOW COUPLER Device and System have been specifically designed for use in end-to-end anastomosis of blood vessels and the detection of blood flow at the anastomotic site. On an as needed basis, blood flow can be detected for up to 7 days.

The FLOW COUPLER System consists of a FLOW COUPLER Device and a FLOW COUPLER Monitor. The FLOW COUPLER Device includes a 20MHz ultrasonic Doppler transducer (probe) attached to one of the FLOW COUPLER rings, and an external lead. The FLOW COUPLER rings are made of high density polyethylene and surgical grade stainless steel pins. A protective cover and jaw assembly protect the rings and probe which allow for easy loading onto the Anastomotic Instrument. Both the protective cover and jaw assembly are disposable.

Accessories to the FLOW COUPLER System include a reusable Anastomotic Instrument (surgical-grade stainless steel and titanium), a reusable Vessel Measuring Gauge (surgical-grade stainless steel), COUPLER Forceps (surgical-grade stainless steel), and a Sterilization Tray (anodized aluminum).

#### INDICATIONS FOR USE:

The FLOW COUPLER Device is a single use, implantable device that is intended to be used in the end-to-end anastomosis of veins and arteries normally encountered in microsurgical and vascular reconstructive procedures. The FLOW COUPLER Device includes a pair of permanently implanted rings which secure the anastomosis and a removable Doppler probe that is press-fit onto one of the rings. When the FLOW COUPLER Device is used in conjunction with the FLOW COUPLER Monitor, the FLOW COUPLER System is intended to detect blood flow and confirm vessel patency intra-operatively and post-operatively at the anastomotic site. Post-operatively, blood flow can be detected on an as needed basis for up to 7 days. The FLOW COUPLER Doppler probe is not intended to be a permanent implant and should be removed 3 to 14 days post-operatively.

#### CONTRAINDICATIONS:

The FLOW COUPLER Device is not indicated for use in end-to-side anastomosis or for patients presenting conditions that would normally preclude microvascular repair with suture technique. Examples of such conditions include, but are not limited to:

- Pre-existing or suspected peripheral vascular disease,
- Ongoing irradiation of the area of reconstruction,
- Clinical infection of the area of reconstruction,
- Anticipated infection due to significant contamination of the area of reconstruction,
- Friability of the vascular tissue due to sclerotic conditions,
- Concurrent diabetes mellitus, or
- Concurrent corticosteroid therapy

The FLOW COUPLER Device and System is contraindicated for use in the central circulatory system.

#### WARNINGS:

- Failure to use the Vessel Measuring Gauge to approximate the vessel size could result in using a FLOW COUPLER Device of an inappropriate size. Using a ring too large for the vessel may result in stressing or tearing of the vessel wall and a compromised anastomosis. Using a ring too small for the vessel may unduly constrict the vessel and lead to thrombosis or ring separation.
- Failure to squeeze the FLOW COUPLER jaws with a hemostat or similar instrument prior to ejection of the joined rings may result in an inadequate friction fit and possible ring separation. **Inspect the anastomotic site** to ensure that the anastomosis has been satisfactorily completed.
- The FLOW COUPLER Device is supplied sterile and is single use only. **Do not resterilize or reuse** the FLOW COUPLER Device.

- Resterilization may compromise the structural integrity of the product which may lead to incomplete anastomosis.
- Device cannot be reused due to possible structural damage during first use which may lead to incomplete anastomosis.
- Do not use the FLOW COUPLER Device if the package has been opened or appears to be damaged or compromised as sterility may be compromised. Failure to observe this warning may result in surgical infection.
- Safe use of the FLOW COUPLER Device for the anastomosis of tubular structures other than veins and arteries has not been established.
- Safe use of the FLOW COUPLER Device for the anastomosis of growing vessels in children or adolescents has not been established. Not intended for fetal use.
- Safe use of the probe portion of the FLOW COUPLER Device during MRI procedures has not been established. Therefore the probe should be removed prior to a MRI procedure.
- Security of an anastomosis utilizing FLOW COUPLER Devices that have been approximated, reopened, and then reapproximated has not been demonstrated. When reapproximation of the anastomosis is desired, the vessel should be removed from each ring and a new FLOW COUPLER Device utilized.
- Ensure that suture sleeve and connectors are not implanted.
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray **must be sterilized prior to use.**
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray should be thoroughly inspected before use. Instruments that are damaged and/or in need of repair should not be used.

#### CAUTIONS:

- Use of the FLOW COUPLER Device involves potential risks normally associated with any implanted device, e.g., infection, perforation, or laceration of vessels, erosion, implant rejection, or device dislodgement/migration.
- The angle of the probe wire relative to the flap will be influenced by the orientation of the Anastomotic Instrument during formation of the anastomosis. To avoid unwanted kinking or twisting of the vessel during positioning of the flap-which may result in poor flap perfusion-care should be taken to establish the desired angle of the probe wire relative to the flap and to adjust the Anastomotic Instrument accordingly prior to starting the anastomosis.
- Should a probe be prematurely removed from the probe-holder, do not attempt to re-insert the probe into the probe-holder. Instead remove rings and implant a new FLOW COUPLER Device.
- **Probe wire is delicate. The use of crushing forceps may damage the probe wire.**

- **Use caution when manipulating the probe wire. Sharp bends may cause damage to the probe wire.**
- **The use of clamps on the external lead wire may damage the external lead.**
- The probe is not intended to be a permanent implant and should be removed 3 to 14 days post-operatively.
- Avoid excessive force to remove the probe from the patient, which may cause injury to the blood vessel. If the probe can not be removed using gentle traction, the probe should be surgically removed. **Do not cut the probe wire.**
- Assure that the probe is attached to the probe wire upon removal of the probe. If not, surgical removal of the probe is required.
- The FLOW COUPLER Device should only be used with the GEM FLOW COUPLER Monitor.
- During the use of all ultrasound devices, the operator should minimize the exposure of ultrasound energy to the patient using the principle of ALARA (As Low As Reasonably Achievable).

#### INSTRUCTIONS FOR USE:

These *Instructions for Use* are designed for proper use of this device. They are not intended to serve as a reference to surgical technique, to supersede institutional protocols or professional clinical judgment regarding patient care.

It is the responsibility of the clinician to inform the patient that he/she is the recipient of permanent implants which contain metal components (surgical-grade stainless steel pins). The FLOW COUPLER rings have been evaluated with a 1.5 Tesla and 3.0 Tesla magnetic field and no change in displacement was observed in each of three orthogonal planes.<sup>1,3</sup> The stainless steel pins in the FLOW COUPLER Devices are nominally nonferromagnetic. However, the US Food and Drug Administration (FDA) has made recommendations for any medical device implanted which have metallic components to include:

- Documentation in the official medical record of the identity of the implant (manufacturer, model number, lot and serial numbers, and identifying marks, if any).
- Documentation of the technique and results of any magnetic testing performed on the implant or that no such testing was done.
- Patient education regarding the particular implant and recommendation for identifying medical alert card, bracelet, or necklace characterizing the implanted device.<sup>2</sup>

#### 3.0MM COUPLER SIZE OR SMALLER:

##### END-TO-END ANASTOMOSIS:

Using conventional microsurgical technique, mobilize a minimum of 1 cm of each vessel end. Using vascular clamps, clamp off the vessel(s) and irrigate the vessel openings. The

FLOW COUPLER Device requires a greater amount of free vessel within the clamps than a conventional suture repair.

1. After gentle dilation, estimate the **outer** diameter of each vessel using the Vessel Measuring Gauge. The circular guides on the gauge **should not** be placed inside the vessel lumen (See Figure 1). If there is a size discrepancy between the two vessels, use the measurement of the smaller vessel to choose the appropriate FLOW COUPLER Device. The degree of vessel spasm and the elasticity of the vessel should be considered when choosing the FLOW COUPLER Device size to be used.
  2. Select the appropriate size FLOW COUPLER Device. Both vessel ends should be approximately the same size as the inside diameter of the FLOW COUPLER Device being selected.
  3. Remove the lid from the outer tray and aseptically remove the inner tray; the inner tray may be placed in the sterile field. Inspect the inner tray. Do not use if the inner tray is damaged or if the seals are not intact. Remove the lid from the inner tray.
  4. Turn the Anastomotic Instrument knob fully counterclockwise, and then insert the FLOW COUPLER Device onto the Anastomotic Instrument while FLOW COUPLER Device is still in tray. **The matching indicator arrows on the FLOW COUPLER Device and the Anastomotic Instrument should be pointing toward each other when loading** (See Figures 2 & 3). **Ensure that an audible click is heard for proper loading.**
  5. Remove FLOW COUPLER Device from tray and protective cover, being careful not to pull the wire (See Figure 4).
  6. Verify probe function by connecting probe to Monitor and irrigating attached probe tip with sterile saline. (Refer to the Flow Detection section of these Instructions for Use for proper connection instructions.) An audible signal from the Monitor verifies proper function of the device. If no signal is identified, refer to the Troubleshooting section of these Instructions for Use.
  7. Visually inspect to see that both rings are seated at the bottom of the U portion of the jaw (See Figures 5a & 5b) and the pins are not bent. If pins are bent, do not attempt to straighten. Instead use a new FLOW COUPLER Device.
- NOTE: To avoid unwanted kinking or twisting of the vessel during positioning of the flap, care should be taken to establish the desired angle of the probe wire relative to the flap and to adjust the Anastomotic Instrument accordingly prior to starting the anastomosis.
8. Place the Anastomotic Instrument perpendicular to the vessel(s), with the FLOW COUPLER jaw assembly near the two vessel ends. Pull one vessel end through one of the FLOW COUPLER rings using microsurgical forceps (See Figure 6). Care should be taken to avoid twisting of the vessel.
  9. Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale onto one pin. Proceeding in a triangular fashion, impale the vessel firmly upon every other pin, completing three pins (See Figure 7). Complete vessel

placement on the ring by impaling the vessel upon the remaining three intermediate pins (See Figure 8). Ensure that both the vessel wall and the intimal layer are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel (See Figure 9).

10. Repeat Steps 8 and 9 to impale the other vessel end upon the second FLOW COUPLER ring.
11. When both vessel ends have been suitably impaled, visually inspect to ensure that both rings are seated at the bottom of the U portion of the jaw (See Figures 5a & 5b) and the pins are not bent. Bring the rings together (See Figures 10 & 11) by turning the Anastomotic Instrument knob clockwise.
12. **Prior to ejecting the joined rings, gently squeeze the end of the apposed jaws with a small hemostat** (See Figure 12) **to ensure ring approximation and a tight friction fit.** Turn the Anastomotic Instrument knob further clockwise to eject the joined rings.
13. Check the anastomosis under the operating microscope before opening the vascular clamps. Remove the clamps and **inspect the anastomotic site to ensure that the anastomosis has been satisfactorily completed** (patent vessel without leakage).
14. To remove the jaw assembly turn the Anastomotic Instrument knob fully counterclockwise (See Figure 13). Press the release button, located near the arrow on the Anastomotic Instrument, and remove the jaw assembly (See Figure 14).
15. Rinse the Anastomotic Instrument with water after use.

### **3.5MM COUPLER SIZE OR LARGER: END-TO-END ANASTOMOSIS:**

1. to 8. Follow the same directions as for 3.0mm FLOW COUPLER Device Size or Smaller End-to-End Anastomosis (Steps 1 through 8).
9. Take a bite of approximately one to two pin diameters of the vessel wall and intimal lining, evert 90 degrees and impale onto the pin situated nearest to the open part of the Jaw Assembly (open end of the U portion of the jaw). Impale the opposite side of the vessel opening to the pin directly across from the initial pin. Next, impale the vessel onto the pins located near the sides of the ring, keeping the vessel as evenly spaced as possible between the four pins (See Figure 15). Continue vessel placement on the ring by impaling the vessel onto the two remaining pins near the open end of the Jaw Assembly. Complete by impaling the vessel onto the last two pins near the bottom of the Jaw Assembly (bottom of the U portion of the jaw); this final step prevents the ring from sliding out of the Jaw Assembly prematurely (See Figure 16). Ensure that both the vessel wall and the intimal lining are fully impaled upon each pin to reduce the risk of thrombosis. Should the vessel wall tear during impalement, remove the vessel, trim the end, and repeat the procedure. For examples of improper impalement of the vessel see Figure 17.

10. Repeat Step 9 to impale the other vessel end upon the second FLOW COUPLER ring.
- 11 to 15. Follow the same directions as for 3.0mm FLOW COUPLER Device Size or Smaller End-to-End Anastomosis (Steps 11 through 15).

### **FLOW DETECTION:**

Prior to closure of the surgical site verify detection of blood flow.

1. Temporarily secure the probe wire to the skin to prevent the weight of the metal connectors from pulling on the probe.
2. Join the probe connector to either end of the external lead. Attach the other end of the external lead to the FLOW COUPLER Monitor.
3. Turn on the FLOW COUPLER Monitor.

NOTE: The FLOW COUPLER Monitor (GEM1020M-2) can be powered by a rechargeable internal Lithium ion battery pack or with the external power supply. If the battery level is low or critically low, use power supply.

The FLOW COUPLER Monitor (GEM1020M) can be powered by batteries (8 AA) or with the external power supply. If the low battery light illuminates, either replace all 8 batteries or use power supply.

- NOTE: For further instructions, refer to the GEM FLOW COUPLER Monitor Instructions for Use.
4. Select appropriate channel on FLOW COUPLER Monitor and listen for blood flow. Adjust volume as needed. If a strong audible signal is not identified, irrigate the site where the probe tip meets the vessel with saline. During irrigation, an audible signal from the monitor verifies proper function of the device.

NOTE: Do not attempt to adjust probe location.

5. When routing wire away from the anastomotic site, a loose suture may be placed around the wire to ensure that it does not affect the orientation of the joined FLOW COUPLER rings. Optimal wire position would be aligned with probe tip (See Figure 18). Do not bend probe wire at a sharp angle. (see Figure 19) See Figure 20 for an example of proper probe wire angle. Carefully position the probe wire to leave enough wire length in the wound, providing slack to assure there is no tension on the anastomosis.
6. Once satisfied with wire placement, use a tack suture on the probe wire at the wound margin (5-0 or similar). Secure the suture sleeve to the skin (suture, tape or staple). Ensure adequate slack in the wire.
7. Following verification of probe function and wire placement, close the incision using standard techniques. Cover exposed probe wire with medical dressing.
8. On an as needed basis, blood flow can be detected for up to 7 days. The probe is not intended to be a permanent implant and should be removed 3 to 14 days post-operatively.

9. When monitor is not being used to detect flow, external lead may be disconnected from the probe by pulling probe connectors apart.

NOTE: Ischemia or reperfusion rate may delay or affect the initial Doppler signal.

NOTE: If blood flow is not detected with the Monitor post-operatively, rely on clinical indications for patient status.

NOTE: Doppler signal may vary during monitoring period.

10. To remove the probe, first detach the suture sleeve and wire from the skin (remove suture, tape or staple). Remove the probe by applying gentle traction to the wire while applying counter pressure externally at the site of incision until the probe is extracted. Inspect to ensure that probe tip is fully intact. If probe is not present, surgical removal is required.

Devices: Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray.

## SPECIAL INSTRUCTIONS

### WARNINGS

- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray are supplied non-sterile and must be sterilized prior to use
- The Anastomotic Instrument, Vessel Measuring Gauge, COUPLER Forceps, and Sterilization Tray should be thoroughly inspected before use. Instruments that are damaged and/or in need of repair should not be used.

### LIMITATIONS ON REPROCESSING

No particular limitations

INSTRUCTIONS		
Point of use:	Rinse all instruments with water after use.	
Preparation for cleaning:	<ol style="list-style-type: none"> <li>1. Using a neutral (pH 7-10) detergent, wash each tool clean of all blood and debris after every use. Scrub each tool with a soft brush. Pay particular attention to areas where debris can accumulate.</li> <li>CAUTION: Use of a cleaner with a pH greater than 10 will remove the anodized layer of the Anastomotic Instrument and the Sterilization Tray.</li> <li>2. Avoid use of any harsh material that can scratch or mar the surface of the instruments</li> <li>3. Rinse the instruments thoroughly with running water. Apply a fine jet stream through the hole in the Anastomotic Instrument knob end and press the release button while rinsing to ensure that all surfaces of the instrument are cleaned.</li> </ol>	
Cleaning: Automated	Using an automated washing machine, clean at a temperature of 45° – 55° C, with a neutral (pH 7-10) cleaning solution for at least 10 minutes.	
Cleaning: Manual	Place the instruments in an ultrasonic cleaner utilizing a neutral (pH 7-10) cleaning solution and clean ultrasonically for 15 minutes. Rinse the Anastomotic Instrument thoroughly again, applying a fine jet stream of water through the hole in the Anastomotic Instrument knob end.	
Disinfection:	(Optional) Using an automated washing machine, thermally disinfect the Instrument at a temperature of 90° - 95° C, for a minimum of 5 minutes.	
Drying:	Following cleaning by either manual or automated cleaning methods ensure instruments are fully dry. Do not exceed 100° C for 30 minutes.	
Maintenance, Inspection and Testing:	<ul style="list-style-type: none"> <li>• Ensure that all visible debris is removed to assure the continued quality of the instruments.</li> <li>• Lubricate the cleaned Anastomotic Instrument (including knob) with a water-soluble lubricant prior to sterilization. Failure to clean and lubricate the Anastomotic Instrument as directed may result in instrument failure.</li> </ul>	
Packaging:	Package the instruments using the appropriate method for the sterilization cycle chosen.	
Sterilization:	<b>GRAVITY AUTOCLAVES</b>	
	Temperature	Recommended Exposure Time (not Total Cycle Time)
	270°F (132°C)	15 minutes (wrapped)
	<b>PREVAC AUTOCLAVES</b>	
	Temperature	Recommended Exposure Time (not Total Cycle Time)
	270°F (132°C)	4 minutes (wrapped)
273°F (134°C)	3 minutes (wrapped)	
	NOTE: It is recommended that each institution establish the efficacy of its sterilization procedure.	
Storage	Recommended storage at controlled room temperature 20-25°C (68-77°F).	

EN

## PROBE TROUBLESHOOTING GUIDE

SYMPTOMS	POSSIBLE PROBLEMS	SOLUTION
No sound output Intra-operatively	No moisture contact	Site irrigation
		Verify no vessel stricture
		Check blood flow (vein and artery)
		Massage blood vessel to increase blood pressure
		Wait until blood flow can be seen and heard
No power	No power	Ischemia or reperfusion rate may delay or affect the initial Doppler signal. Check with hand held probe.
		Check all connections: • Probe Connector to External Lead • External Lead to Monitor
		Replace External Lead
		Rely on clinical indications for patient status.
No sound output Post-operatively	No power	External Lead is not functioning
		Probe is not functioning
		Check all connections: • Probe Connector to External Lead • External Lead to Monitor
		Rely on clinical indications for patient status.
No sound output Post-operatively	External Lead is not functioning	Replace External Lead
		Probe may have lost contact with vessel
		Rely on clinical indications for patient status.
No sound output Post-operatively	Probe is not functioning	Rely on clinical indications for patient status.

## MONITOR TROUBLESHOOTING GUIDE

SYMPTOMS	POSSIBLE PROBLEM	SOLUTION
No sound output	No power	Verify Monitor power is on
	No power	Check all connections: • Probe Connector to External Lead • External Lead to Monitor • Monitor to AC Power Supply • AC Power Supply to Power Cord • Power Cord to outlet
	Volume is too low	GEM1020M: Adjust volume using Volume Increase switch GEM1020M-2: Adjust volume using Volume Increase Button or LCD volume control
	Wrong channel is being used	GEM1020M: Verify the correct channel is illuminated GEM1020M-2: Verify the correct channel button is illuminated and the correct channel is displayed on the LCD screen
	GEM1020M: Batteries are dead GEM1020M-2: Critically Low Battery screen displays on LCD screen	GEM1020M: Replace batteries or use external power supply GEM1020M-2: Connect the Power Supply to the monitor
	Monitor not functioning	Connect a different Monitor
	Monitor not functioning	Contact Synovis Micro Companies Alliance
	Weak sound output	(GEM1020M only): Weak batteries (low battery indicator is illuminated)
Volume is too low		Adjust volume using Volume Increase switch
Monitor not functioning		Connect a different Monitor
Monitor not functioning		Contact Synovis Micro Companies Alliance
(GEM1020M-2 only): LCD touch screen malfunction	LCD touch screen unresponsive	Power off and on the monitor Contact Synovis Micro Companies Alliance Customer Service

**SERVICE:**

For Customer or Technical service, contact:

Phone: + 205.941.0111 or 1.800.510.3318

Fax: + 205.941.1522

Website: synovismicro.com

**ACCESSORIES & PARTS:**

ITEM	REF
FLOW COUPLER Monitor	GEM1020M, GEM1020M-2
Power Supply	GEM 1020PS, GEM1020PS-2
FLOW COUPLER Device 2.0	GEM2752-FC
FLOW COUPLER Device 2.5	GEM2753-FC
FLOW COUPLER Device 3.0	GEM2754-FC
FLOW COUPLER Device 3.5	GEM2755-FC
FLOW COUPLER Device 4.0	GEM2756-FC
FLOW COUPLER External Lead	GEM1003EXT-FC

**DISCLAIMER OF WARRANTY:**

Synovis Micro Companies Alliance, Inc., a subsidiary of Baxter International Inc. warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. Since SMCA has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SMCA does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SMCA will replace any device which is defective at the time of shipment. No representative of SMCA may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

**REFERENCES:**

1. DeLacure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

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#### LITERATUR:

1. DeLaCure M and Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
2. Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Number 2, June 1993.
3. Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

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#### ERWIJZINGEN:

- DeLacure M en Wang H: Magnetic Resonance Imaging Assessment of a Microvascular Anastomotic Device for Ferromagnetism. Journal of Reconstructive Microsurgery 13:8, 1997.
- Caution needed when performing MRI scans on patients with aneurysm clips. FDA Medical Bulletin Volume 23, Nummer 2, Juni 1993.
- Synovis MCA Internal report, V-0718R rev A, Coupler MR Compatibility Verification Report.

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Flow Coupler and the GEM design are trademarks of Baxter International Inc., or its subsidiaries.
Synovis Micro Companies Alliance, Inc. (a subsidiary of Baxter International Inc.),