World of Microsurgery  
*Welt der Mikrochirurgie*

CO

**S&T**

# Contents

## *Inhalt*

---

**World of Microsurgery** 0.01-06  
*Welt der Mikrochirurgie*

---

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Index  
*Index*

# S & T<sup>®</sup>

All **S&T** products must be used exclusively for their intended purpose.

All instruments are subject to technical changes in design and construction. We will supply the current model.

Our general terms and conditions of sale and deliveries are applicable.

Reprints of this catalogue in whole or in part is not permitted. Printed or electronic matters of any kind for which illustrations of our products are required with or without the respective article number are subject to our written consent.

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History  
*Geschichte*

**0.02**

Thanks to the Surgeons  
*Dank an die Chirurgen*

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## To our Customers *An unsere Kunden*

### Dear Customer...

welcome to S&T's latest catalog — your resource and authority in the field of Microsurgery.

#### Microsurgery

...literally the cutting edge of the surgeon's art, has been practiced in certain procedures for well over 50 years.

Nonetheless, it became truly viable with the advent of microsurgical needlesuture combinations which, in the smaller dimensions, were first manufactured by S&T in the late 60's of the last century.

Working together with world renowned surgeons, Werner Spingler and Eugen Tritt, the founders of S&T, manufactured even better and finer needles which are used all over the world by surgeons, to the benefit of their patients.

However, it was soon noticed that an essential ingredient was missing, and surgeons started asking for modified instruments, better capable of handling the fine needles and structures.

In 1980, S&T's involvement in instruments had increased to such proportions that two spin-off companies were founded to efficiently handle these requests.

Since then, S&T has undergone various transformations, which culminated in the merger of the S&T companies in the S&T AG in the new millenium.

As of July 2021 S&T is a subsidiary of the Swedish Stille Group and proud to now have the opportunity to develop and grow together in both premium segments.

Today, S&T is directly represented in most major microsurgical markets with state-of-the-art product lines ranging from forceps, to needle holders, scissors, vessel clamps as well as special instrumentation.

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**S&T**

# To our Customers *An unsere Kunden*

S&T is proud to acknowledge its collaboration with the world's most prestigious microsurgions. Together, we have crafted some of the most innovative products in the field. S&T products are crafted from premium materials, and are built from a sound foundation of research, expert cooperation, and precise engineering. We at S&T are so confident in our products and processes that we challenge you to find a better line of instruments anywhere else in the world.

**All S&T instruments carry a lifetime warranty against defects in workmanship and materials.**

We thank you for your loyalty and support. Our pledge to you is that S&T will continue its dedication to excellence — so you can achieve it.

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## To our Customers *An unsere Kunden*

### Service partners

S&T has service partners all over the world. You may contact our partners for answers to questions relating to our products and services (such as workshop and training course support for microsurgical technique, or teaching/training materials for the care and handling of surgical instruments). For a list of our business partners and products as well as services we refer you to our web site

[www.microsurgery.ch](http://www.microsurgery.ch)

Of course you may also contact us directly.

**Phone: +41 52 672 74 07**  
**Fax: +41 52 672 74 02**  
**E-Mail: [info@s-and-t.net](mailto:info@s-and-t.net)**  
**S&T AG**  
**Tobelraastrasse 2**  
**CH-8212 Neuhausen**  
**Switzerland**

**E-Mail: [europe@s-and-t.net](mailto:europe@s-and-t.net)**  
**S&T Medical GmbH**  
**Bachstrasse 3**  
**D-79798 Jestetten**  
**Germany**

It will be our pleasure to see to all your needs in instrumentation. Modifications and ideas for new instruments are always welcome.

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# To our Customers An unsere Kunden

## Our contribution to quality

Since 1966 S&T is developing and manufacturing outstanding microsurgical sutures and instruments for mirosurgeons around the globe. Technical capability in combination with innovation, reliability and a sophisticated measurement and monitoring system are the cornerstones of our integrated approach focused on quality. To manufacture a state of the art product we take worldwide variably surgical techniques as well as all special requirements of the end-user into consideration. Together with them we constantly develop and manufacture customized solutions bearing all relevant standards in mind.

Together and in partnership with your Notified Body TÜV Product Service in Munich/Germany we are able to confirm the implementation of all relevant national and international standards. The fulfillment of all regulatory requirements of our QM system according to EN ISO 13485: 2016 is a matter of course as well as the compliance with relevant standards regarding packaging, sterilization, risk management and clinical trials.

All S&T products undergo a **100% final inspection** to avoid even the smallest statistical breakout. Therefore we take pride in issuing a life time warranty on material and craftsmanship on all our instruments.

In this spirit we constantly work on giving you as an enthusiastic and loyal S&T customer the pure conscience to recommend us to others.

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

## Care of Instruments

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## Instrument Care and Handling

Proper care and handling techniques will help preserve the life of an instrument, and ensure that it functions at peak performance. Avoid having your instruments lay loose together. You should not combine delicate micro instruments with macro instruments as this may cause damage to the delicate tips. Most damage to microsurgical instruments occurs during processing and cleaning, so it is crucial that you handle the instruments with great care. Finally, we recommend you invest in a sturdy instrument container that is specifically designed for the storage, sterilization and transport of microsurgical instruments.


**Important**  
Most instrument damage occurs during processing and cleaning. The best way to avoid this is by handling them with great care. This is especially true for microsurgical instruments.

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### Pre-cleaning

It is best to remove excess tissue or other residues from your instruments as soon as possible, ideally immediately after the procedure. This can be done manually or by an ultrasonic cleaner. We suggest that a disinfecting agent be used. Keep in mind that temperatures below 40°C are needed to avoid protein coagulation.

Make sure to move all approximators, ratchets and sliding parts to ensure that all residues are removed. Clamps, scissors and needleholders need to be opened to allow contact with cleaning solutions.

### Cleaning

In general, you should follow the pre-scribed guidelines established for your institution. It is recommended that cleaning and disinfecting be done in a thermal disinfectant. Initial cleaning commences at temperatures of approximately 40-60°C with pH-neutral or alkaline cleaners. Rinsing is best done with cold water and without additives. Foam generation should be avoided as it may inhibit the cleaning action.

You should follow the manufacturers' instructions regarding use of cleaning agents and cleaning equipment. Water used in all cleaning processes should have a low concentration of chloride ions, as excessive concentrations may lead to pitting and other forms of corrosion. Final rinsing should be done with distilled / deionized water. This helps to avoid staining of the instruments.

#### Important

**To avoid damaging the delicate tips of micro instruments, process them in their own designated tray. Avoid overloading the tray, making sure that micro instruments do not come in contact with macro instruments. Scissors and needleholders should be opened to ensure proper cleaning.**

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# S&T

## Care of Instruments

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Be mindful of instruments with large surface areas. They can overshadow smaller instruments situated below them, inhibiting them from being cleaned properly and effectively.

Instruments should be dried thoroughly after cleaning, to avoid corrosion.

Only clean instruments should be sterilized.

#### **Important**

**New instruments should be processed and sterilized prior to their first use.**

### **Disinfection**

For the disinfecting process, we recommend using a cycle with the highest temperatures (usually 80-95°C).

### **Inspection**

Instruments should be inspected after each use and immediately after cleaning. The instruments tips should approximate well. Ensure that all scissors and needle-holders are in proper working order. Pay special attention to hinges and box locks. Make sure that they are thoroughly cleaned and function smoothly. Check for bent and broken tips. Instruments that are damaged should be removed from service and either be repaired or replaced.

Instruments that are beyond repair should be discarded, ensuring that they will not find their way back into the operating room.

If corrosion is apparent, the instrument should be removed from service and sent out for repair. Depending on the nature of the corrosion, some mild forms may be polished away. In the case of serious pitting or corrosion, the instrument should be discarded as to avoid cross contamination with other instruments.

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### Lubrication

#### **Instruments with moving parts perform best when properly lubricated.**

Instrument lubrication not only makes an instrument perform better, but also reduces the effects of corrosion, ultimately prolonging the life of the instrument. We recommend any paraffin-based lubricant that is certified for use with medical devices. Remember to use only enough to wet the relevant areas, allowing instruments to drip dry before placing them in a tray or wrapping.

S&T instruments have a **lifetime warranty** on workmanship and materials. Instruments that are judged to be defective should be returned to S&T. Contact your local representative. When returning an instrument, supply as much of the following information as possible: Article code or catalog number; numerical code, located next to the CE mark; and description of the defect. This information is used to improve S&T instruments and ensure that our products remain state-of-the-art.

### Packaging

Packaging plays an important role in the sterility of your instruments. Corrosion may occur if instruments are handled improperly during the packaging process, as the drying cycle in the sterilizer may be inhibited. Further, it is possible to invalidate the entire sterilization process if defective packaging materials do not properly protect the sterilized instruments.

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# Care of Instruments *Instandhaltung der Instrumente*

## Sterilization

In general, you should follow your institution's set guidelines and procedures when sterilizing microsurgical instruments. The most common form of sterilization is steam autoclaving. Microsurgical instruments should be sterilized in their own case. We recommend autoclaving your S&T instruments in saturated steam, at a minimum of 132-134°C, for 3-18 minutes. Stained instruments are an indication of problems having to do with cleaning/sterilization processes relating to one or more of the following:

- cleaning agents
- steam quality in sterilization
- water quality
- other agents

**Important**  
**Non-stainless or corroded instruments may damage other instruments during sterilization.**

## Storage

Storage sites should be clean, dust-free and well ventilated. Environmental parameters such as temperature, humidity and radiation should be congruent with the standards for human habitation. Humidity must be below the dew point, at all times.

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## Care of Instruments

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#### Use

The life of an instrument may be extended if it is cleaned as often as possible during a procedure. Blood and other contaminants should be removed from the surface of the instrument before it is allowed to dry. This is best done with distilled /deionized water, **not** with saline or ringers solution.

Instruments that appear or become damaged during a procedure should be removed from use and marked for later evaluation, repair and/or replacement.

An influencing factor in an instrument's functionality is the transportation of the instrument between processing facilities, storage and use. We recommend the user evaluate methods and establish a system that is as simple and safe as possible.

#### Important

The guidelines that we have provided here are general in nature. Care and handling procedures vary from country to country and from institution to institution. You should adhere to procedures established by your local protocol. If you require additional information please visit our homepage [www.microsurgery.ch/downloads](http://www.microsurgery.ch/downloads) for a pdf-download of our validated cleaning processes.

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
### Grund-Information

## Scissors

All S&T microsurgical scissors are skillfully handcrafted to microsurgical scissors feature: exacting standards. You will immediately notice the advantages of these scissors:

- Flawless design due to years of experience in the field
- A uniquely **smooth cutting action** due to precise engineering
- **Minimum closing pressure** in operation due to equal pressure distribution along the blades

## S&T® Microsurgical Scissors Patterns

### Dissection

This pattern features **curved** blades with tips that are rounded to avoid possible damage to structures during dissection.

**SDC = scissors, dissection, curved**

### Adventitia

This pattern features **straight** blades with tips that are pointed for fine trimming and suture cutting.

**SAS = scissors, adventitia, straight**

### Nerve/Tendon

This pattern features straight blades with one **finely serrated** blade, ensuring a controlled even cut when preparing tissue such as nerves or tendons. Due to design characteristics, S&T serrated scissors provide exceptionally smooth cutting action compared to other serrated scissors.

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### Grund-Information

## Dissecting or Adventitia Pattern?

Microsurgical scissors in the S&T product offering are available in either dissecting or adventitia patterns. The distinction lies in the sharpness of the tips and the curvature of the blades. Dissection scissors are rounded at the tip with a radius of 0.3 mm, with blades that are slightly curved on the flat. Adventitia scissors are very sharp at the tip, with straight blades.

Dissection scissors are used for separating nerves and blood vessels from the surrounding tissue, without damaging them. When held closed, they are often also used as a dissecting probe. Dissection scissors are unsuitable for fine trimming work, since at high magnification the rounded tips prevent the surgeon from seeing exactly the tissue that is being cut, at the moment of cutting.

Adventitia scissors are used mostly at high magnification, for fine trimming and re-shaping of structures that have already been dissected free. They are also used for cutting fine micro suture material, 8/0 or smaller.

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## SDA-type scissors

SDA-type scissors were specifically developed for fine dissection of the extremities and similar applications. Unlike other dissection scissors, the blades are **straight**. These scissors are held between thumb (top position) and index, stabilized with the middle and ring finger. You will immediately notice that handling the scissors in this position is easier than the usual "chuck grip". A further advantage is the better view of operating field.

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## Basic informations Grund- Information

### Blade lengths

The blade length is measured from screw to the tip.

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### Westcott scissors

S&T Westcott scissors feature finer blades than regular Westcotts, originally used in plastic surgery. The typical shape of the blades remains.

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### Serrated scissors

The unique design of these handcrafted S&T scissors ensures very smooth cutting action, even with the serrations.

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### Vannas scissors

Originally designed for ophthalmic procedures, the vannas scissors are applied in vascular surgery, too. They feature finer, smaller and shorter blades.

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


**Scissors by Blondeel**  
*Scheren nach Blondeel*

**Round handle**  
*Rundgriff*

These scissors were designed by Professor Phillip Blondeel for use in breast reconstruction utilizing a perforator flap. While the ring handle configuration makes these suitable as preparation scissors for fine structures, the round handle, spring instrument configuration, with the curved blades, is extremely useful for fine trimming and dissection work. This instrument will perform tasks for which you might normally need different instruments.

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


Straight Dissection or curved Adventitia scissors? Our specially designed scissors meet all specific requirements!

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



Index  
*Index*

**S&T**  
**Basic**  
**informations**  
*Grund-*  
*Information*

All S&T forceps are handcrafted under the microscope, to the highest standards. Each S&T forceps is inspected utilizing rigid QC standards, assuring you consistently well made instruments, designed specifically for microsurgery.

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Flat handle  
 Classic line  
 Flachgriff  
 Classic line

Round handle  
 Balanced line  
 Rundgriff  
 Balanced line

**Advantages of Balanced Instruments**

Balanced instruments are designed to put the center of gravity (the balance point) in the web space between the thumb and index finger. This is accomplished by adding a slight counter weight to the proximal end of the instrument. The design of balanced instruments enhances tactile feel, aiding in better control and precision. Balanced instruments are less tiring to handle, especially during lengthy procedures. Round handle balanced instruments have better handling characteristics than the standard (non-balanced) round handle instruments.

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**Platform Forceps**

Platform ("Plateau") forceps were initially developed for the US market. The plateau depicts the usable area with "guaranteed" parallel closure. In this way a reliable grip on sutures can be ensured.

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
**S&T**  
**Basic**  
**informations**  
*Grund-*  
*Information*

Note the extreme adherence to mathematically exact curvatures. The finishing of all S&T forceps under the microscope and rigid QC measures assure extremely well made instruments suited for use even in the most adverse situations.

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
**Unique gripping grooves**

Looking at the details of the gripping grooves it will be noticed that at the top, where groove meets "median" surface, the angle is only slightly beveled to enhance gripping power while ensuring that surgical gloves will not be damaged. However, the bottom of the grooves have a radius, as dirt will have less chance of accumulating and sticking to a radius as opposed to an inside angle. The gradation and depth of the grooves was designed to ensure that the glove material would form "into" the grooves, again maximizing feel on the instrument. The form of the striage has also been chosen because it does not tire the tactile nerves of the hand.

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**Angulated and Curved Forceps**

Many instruments feature angulated or curved tips, giving you a grip that is parallel/perpendicular to the working surface, depending on the object you are handling. While most surgeons are comfortable with our standard 45° or 36° angles, others may prefer the flexibility and curvature of an instrument that conforms to any hand position. Curved forceps are designed to facilitate needle driving or passage through tissue. Straight forceps are commonly held in the non-dominant hand and are used to either grasp tissue or receive the needle when suturing.

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**FRS-15 RM-8 (1:1)**

6" 150 mm	8 mm	0.3 mm	Balanced Line, Plateau	<b>FRS-15 RM-8</b>	00165
6" 150 mm	8 mm	0.3 mm	Balanced Line, Plateau, <b>CRYSTAL</b>	<b>FRS-15 RM-8 TC</b>	00761
6" 150 mm	8 mm	0.1 mm	Balanced Line, Plateau, <b>SUPERFINE</b>	<b>FRS-15 RM-8.1</b>	00762
6" 150 mm	8 mm	0.3 mm	Balanced Line, 1 x 2 teeth <i>Balanced Line, 1 x 2 Zähne</i>	<b>FRSF-15 RM-8</b>	00669
6" 150 mm	8 mm	0.5 mm	Balanced Line, without Plateau <i>Balanced Line, ohne Plateau</i>	<b>FRS-15 RM-8 WP</b>	00936

**FRS-18 RM-8 (1:1)**

7 1/8" 180 mm	8 mm	0.5 mm	Balanced Line, Plateau	<b>FRS-18 RM-8</b>	00235
7 1/8" 180 mm	8 mm	0.1 mm	Balanced Line, Plateau, <b>SUPERFINE</b>	<b>FRS-18 RM-8.1</b>	00361
7 1/8" 180 mm	8 mm	0.3 mm	Balanced Line, 1 x 2 teeth <i>Balanced Line, 1 x 2 Zähne</i>	<b>FRSF-18 RM-8</b>	00670
7 1/8" 180 mm	8 mm	0.5 mm	Balanced Line, without Plateau <i>Balanced Line, ohne Plateau</i>	<b>FRS-18 RM-8 WP</b>	00937
7 1/8" 180 mm	8 mm	0.3 mm	Balanced Line, without Plateau <i>Balanced Line, ohne Plateau</i>	<b>FRS-18 RM-8 WP3</b>	00945

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
**FRC-15 RM-8 (1:1)**

6" <i>150 mm</i>	8 mm	0.3 mm	Balanced Line, Plateau	<b>FRC-15 RM-8</b>	00166
---------------------	------	--------	------------------------	--------------------	-------

**FRC-18 RM-8 (1:1)**

7 1/8" <i>180 mm</i>	8 mm	0.5 mm	Balanced Line, Plateau	<b>FRC-18 RM-8</b>	00238
-------------------------	------	--------	------------------------	--------------------	-------

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**Forceps by De Bakey**  
*Pinzetten nach De Bakey*

**Round handle, straight**  
*Rundgriff, gerade*

The atraumatic forceps by De Bakey are specially designed for a gentle grip of tissue. With a width of the tip of 1.1 mm this forceps is suited for very delicate procedures.

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
Distinguished by their minute "hole in the tip", Pierse forceps were originally designed to hold the sclera — the most slippery tissue in the body. As their remarkable gripping power has been appreciated, their use has gradually spread from eye surgery to wider microsurgical fields. The two minute tissue contact pads, backed by the tissue retention hole, securely grip tissue while minimizing trauma to the tissue. This is much better than the interlocking jaws of a conventional fine-toothed forceps. For grasping and holding tissues as diverse as vascular sheath, epineurium, vas deferens or fallopian tube, Pierse forceps are among the most indispensable instruments of a microsurgeon.

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
**Forceps by Sørensen**

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**Round handle, angulated 45°**

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These forceps can be used in conjunction with multipurpose forceps or alone. These angulated forceps are especially suited for work in the horizontal plane.

**FRAS-15 RM-8 / FRAS-18 RM-8**

Multipurpose forceps by Sørensen:

These designs by Jesper Sørensen, M.D., Ph.D., are useful in all plastic and reconstructive surgery procedures. The working tips are configured to enable their use as dilators at the front but also as needle holders and tying forceps further in the back.

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
This handy instrument, with its delicate curved tip, is designed for end-to-side anastomosis. It is specifically intended to make the work of the surgeon's left hand easier. In passing the needle it provides internal counter-pressure or external edge-eversion, even at those angles where a straight forceps would be of little help. It is also an excellent suture tying instrument.

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These forceps are available in 13.5 and 18 cm flat handle as well as 15 and 18 cm round handle versions.

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These forceps are widely known as "hockey sticks", due to their easily distinguishable form.

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Dr. Joseph Banis, University of Louisville, Ky., USA


*Dr. Joseph Banis, University of Louisville, Ky., USA*

length in inch <i>Länge in mm</i>	width handle <i>Griffbreite</i>	Ø-Tip <i>Ø-Spitze</i>	Remarks <i>Bemerkungen</i>	Art. Code <i>Art. Bez.</i>	Art. No. <i>Art. Nr.</i>
--------------------------------------	------------------------------------	--------------------------	-------------------------------	-------------------------------	-----------------------------

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

These forceps are used to introduce micro catheters into blood vessels for medication and diagnostic purposes.

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They are popular in both clinical and lab environments.

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Tissue Forceps  
*Gewebepinzetten*


Round handle, straight  
*Rundgriff, gerade*

Tissue Forceps with ring 1 mm dia., hole 0.5 mm dia., TC coated.

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TFS-18 RM-8 TC (1:1)

7 1/8" 180 mm	8 mm	1.0 mm	Balanced Line, <b>CRYS&amp;TAL</b>	TFS-18 RM-8 TC	00576
------------------	------	--------	------------------------------------	----------------	-------

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


"The special Vascular Dissecting Forceps of S&T are used as fine dissectors for precision work around the small vessels such as perforators during operation involving perforator flaps and other flaps. These dissectors allow for precision dissection with minimal trauma to the adjacent tissue."

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



Basic  
informations  
*Grund-  
Information*

## Vessel Dilators

S&T Vessel dilators feature highly polished parallel tips. They are designed for controlled, interluminal vessel dilation. S&T offers a wide range of styles to meet your specific requirements.

The special ball point form of the tips ensures that the dilating pressure is evenly distributed over the entire area of the vessel. This dilator form was especially designed for use in microsurgery and is **extremely gentle to small vessels**.

S&T's finishes these forceps to the most exacting specifications, ensuring that form, surface and dimensions are exactly within the parameters defined by the developer surgeons.

The mirror-polished surface of the dilator's working area under magnification. **This ensures that unnecessary traumatic damage to the vessel interior is avoided.** The form of these dilators was developed to open the vessel over the entire spasm-affected area in a parallel fashion. The entire dilated area is opened to the same diameter, and the pressure is distributed evenly.

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
**FRS-15 RM-8d.3 (1:1)**

6" <i>150 mm</i>	8 mm	0.3 mm / 4 mm	Balanced Line, Plateau	<b>FRS-15 RM-8d.3</b>	00590
6" <i>150 mm</i>	8 mm	0.1 mm / 4 mm	Balanced Line, Plateau, <b>SUPERFINE</b>	<b>FRS-15 RM-8d.1</b>	00577

**FRS-18 RM-8d.3 (1:1)**

7 1/8" <i>180 mm</i>	8 mm	0.3 mm / 4 mm	Balanced Line, Plateau	<b>FRS-18 RM-8d.3</b>	00594
7 1/8" <i>180 mm</i>	8 mm	0.2 mm / 4 mm	Balanced Line, Plateau	<b>FRS-18 RM-8d.2</b>	00941
7 1/8" <i>180 mm</i>	8 mm	0.1 mm / 4 mm	Balanced Line, Plateau, <b>SUPERFINE</b>	<b>FRS-18 RM-8d.1</b>	00588

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FRC-15 RM-8d.2 (1:1)

6" 150 mm	8 mm	0.2 mm	Balanced Line, <b>SUPERFINE</b>	<b>FRC-15 RM-8d.2</b>	00919
7 1/8" 180 mm	8 mm	0.2 mm	Balanced Line, <b>SUPERFINE</b>	<b>FRC-18 RM-8d.2</b>	00918

#### **CRYSTAL**

TC = Tungsten carbide tips for better grip.

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
Basic  
informations  
*Grund-  
Information*

## Needleholders

There is more to a good needleholder than meets the eye.

This is especially true for S&T's line:

- Secure grip on even the smallest needle and finest suture material is possible due to parallel closure over the entire jaw length.
- Easy tying of the finest suture material is possible due to perfectly designed jaws.
- Closing pressure is precisely calibrated to facilitate holding of the instrument while minimizing hand fatigue.
- Safe suture handling as all sharp corners and rough spots are eliminated.

## Selection of the proper Needleholder

In addition to the observations described in the beginning pages of the catalog, there are four points that you should take into consideration when selecting a needleholder:

- your preference for the needleholder tip, straight or curved
- your preference for locking or non-locking needleholder
- your preference for a hinged or non-hinged needleholder
- Round or flat handle

## Needleholder tip

Curved tip needleholders are the most commonly used instruments, however on occasion some surgeons prefer straight needleholders.

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## Basic informations

### Grund-Information

Needleholders with the designation "C" in the article code are based on the design by Castroviejo and are basically fine, flat handle needleholders. The designation "B" stands for round handle needleholders which were introduced by Barraquer. Both designs originate from ophthalmic instrument designs.

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*Prof. Dr. Ignacio Barraquer y Barraquer, Barcelona, Spanien*

## Needleholders with lock

The lock on a needleholder is mainly intended to be a convenient means for parking the needle and passing it to the surgeon. The lock is not intended to be an aid for holding the needle while suturing.

S&T needleholders with a lock are carefully designed and calibrated. The lock may be disengaged and re-engaged without the inconvenience of letting go of the needle. Engagement and disengagement is marked by an audible click. Since the finger pressure required to release the lock may cause unwanted movement at a critical moment, we recommend that the passing of the needle through tissue should be done with the lock not engaged.

**Important!**  
Refrain from using microsurgical needle holders to hold a macrosurgical needle. Needles over 8 mm length are so thick and the force required to grip and pass them is so high, that even a single use may permanently damage the jaws of a micro needleholder.

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

### Round handle, straight / curved

*Rundgriff, gerade / gebogen*

Curved forceps are extremely useful light-duty needle holders, very widely preferred for lab work. They are also useful for fine clinical work on vessels up to 2 mm in diameter. Above this size a hinged needle holder is recommended, since a firmer grip on the needle is required when working with thicker tissue.

Needleholder forceps have stronger tips than standard microsurgical forceps. Needleholder forceps are produced to the same rigid specifications of our spring-handled needleholders. S&T needleholder forceps are designed to facilitate rotational technique.

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**NHFS-15 RM-8 (1:1)**

6" 150 mm	8 mm	0.6 mm	Balanced Line, straight <i>Balanced Line, gerade</i>	<b>NHFS-15 RM-8</b>	00277
7 1/8" 180 mm	8 mm	0.6 mm	Balanced Line, straight <i>Balanced Line, gerade</i>	<b>NHFS-18 RM-8</b>	00278

**NHF-15 RM-8 (1:1)**

6" 150 mm	8 mm	0.6 mm	Balanced Line, curved <i>Balanced Line, gebogen</i>	<b>NHF-15 RM-8</b>	00167
7 1/8"	8 mm	0.6 mm	Balanced Line, curved	<b>NHF-18 RM-8</b>	00176

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Double action needleholders are especially advantageous in situations where access to an operation site is limited. Double action needleholders feature a reduced handle width compared to a conventional spring handled needleholder.

These instruments are an original **S&T** development.

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


Prof. Dr. Robert D. Acland, University of Louisville,  
Louisville. Ky. USA

### The Original

What started as a laboratory instrument over thirty years ago, has become the most popular micro clamp due to the many advantages this clamp series has for its users:

- light and compact design
- extremely durable
- corrosion resistant metallurgy
- a unique gripping surface
- smooth sliding bar action on all approximator clamps
- clamp profile with parallel sided jaws
- choice of matte or black finish

Offering a range of 5 sizes, the S&T clinical micro clamp series will meet the most exacting needs of small vessel surgery in the diameter range from 0.4 to 5 mm.

#### Clamp Security

All S&T clamps are precisely calibrated to exacting standards. All S&T clamps undergo a rigorous QC testing procedure. S&T clamps are **100% inspected** prior to shipment. S&T's clinical micro clamps are the gentlest and least traumatic of all small vessel clamps.

**Important!**  
Acland® is a registered trademark of  
S&T AG.

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## Basic informations Grund- Information

### Slippage

The flat and non-serrated "Ultra-Grip" surface of these clamps ensures that the pressure is spread evenly over the entire area of the vessel wall. Minimizing over-pressure peaks, which may lead to damage of the vessel interior.

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## Options


### Clamp size?

In choosing the correct clamp size for your application, please refer to the section "Clamp Size and Clamp Pressure" below.

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### "A" or "V" Pattern?

Throughout the series, the "V" pattern is the all-purpose clamp, commonly used on all veins and almost all arteries. The "A" pattern is designed for exceptionally thick-walled and slippery arteries. The slight incurve of the tips of the "A" pattern jaw provides extra security against slippage.

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### Black or Plain Finish?

All clamps are now available in either a black or matte (dark gray) finish. The black finish helps minimize eye fatigue, and also helps avoid burned-out highlights in photography and on-screen. The matte finish features a dark gray, minimally reflective surface. The selection of a black or matte finish clamp is a matter of preference.

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### Grund-Information


#### Clamp Size and Clamp Pressure

When choosing the right size of clamp for your work, clamp pressure is an important consideration. Clamp pressure is determined by the closing force of the clamp (in grams), divided by the area of the vessel wall that is compressed between the jaws (in square millimeters). Thus, the smaller the area, the higher the pressure is, inversely proportional to the vessel size. Each size of clamp has its own different and carefully calibrated closing force. Each clamp size is designed for a specific range of vessel diameters (Fig. 1).

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Fig. 2

#### Leakage

Even on the largest, recommended diameter vessel (Fig. 2), each size of clamp will withstand blood pressure within the vessel up to 200 mm Hg, without leakage.

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Despite the wide range of closing forces, the pressure exerted by each clamp, big or small, remains consistent throughout the series. Each size of clamp exerts a pressure of 5gm/mm<sup>2</sup> when used on the largest vessel in its range, and 15gm/mm<sup>2</sup> when used on the smallest vessel. Even at the higher range, S&T's clinical micro clamps have the gentlest working pressure of any small vessel clamp available.

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The table (Fig. 1) summarizes the relationship between clamp size, recommended vessel diameter and clamp pressure.

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"Vessel diameter" refers to the external diameter of the vessel in its natural state of full dilation.

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
## Basic informations

### Grund-Information

There is a 50% overlap in the recommended vessel diameter range, between each clamp size and the next larger or smaller size. Thus, for every vessel diameter, with the exception of those at the extreme ends of the range, there will always be two suitable clamps. Whenever a choice arises, the smaller clamp should be selected to ensure a gentler, less pressured handling of the vessel.

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### Our new design

#### Advantages

- Marking shows the appropriate position of vessels
- Border lines show the range of the vessel diameters
  - The outer lines designate the maximum vessel diameter (Max)
  - The inner lines designate the minimum vessel diameter (Min)
- The **S&T** logo: The certainty that you are using original Acland<sup>®</sup> microvascular clamps!

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
**S&T**  
Basic  
informations  
*Grund-  
Information*

## Single clamps

Single clamps are used in situations where space is at a premium and you find it difficult to place an approximator clamp. Generally you will use single clamps when you are working on end-to-side anastomosis and you need safe haemostasis on the flap vessel.

Packaging unit: 2 single clamps

## Approximator clamps

Approximator clamps are the standard clamps for end-to-end anastomosis or if the recipient vessel on an end-to-side anastomosis is to be clamped and approximated. The ease with which vessel ends may be positioned as well as the exceptional access to the procedure site has made these clamps top sellers.

Packaging unit: 1 Approximator clamp

## Care of Clamps

Micro clamps are prone to damage from the accumulation of blood clotting in tight places. With proper care and cleaning, this can be avoided. Blood should not be allowed to dry on the clamps during surgery. For cleaning after surgery, we recommend immersion for 30 minutes in hemolytic detergent. When immersing the clamps, hold them open for a moment to wet the inner surfaces. Be sure to slide approximator clamps along the bar, back and forth a few times, to clear the sliding units. After immersion, rinse vigorously with a fine, high-pressure water jet to remove any blood residue from spring interspaces, sliding units and suture cleats.

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Acland® Micro  
Vascular Clamps  
*Acland® Mikro-  
Gefäßklemmen*

**Approximator clamps  
with frame**

These clamps are generally used in labs or in instances where there is no assistance available. The cleats on the frame are used for belying the suture material, ensuring perfect positioning during an anastomosis.

**Important!**

**Do not use a clamp on a vessel below the recommended diameter range. Use of any clamp on a vessel too small will result in unacceptably high pressure and risk of vessel damage.**

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
# Acland® Micro Vascular Clamps *Acland® Mikro- Gefässklemmen*

## **Applying the Clamps**

The three smaller sizes, B-1, B-2, and B-3, are designed to be applied with the Meyer clamp applicator forceps. The forceps are direct acting. When the forceps is squeezed, the clamp opens. When held in the forceps, the clamp can be manipulated firmly in any plane with-out the danger of slipping or falling out. Clamp applicator forceps are available with or without a lock.

RD and HD clamps are designed to be applied by hand. Their broad, corrugated gripping surfaces make this easy.

Prof. Dr. Viktor E. Meyer, University Hospital Zurich, Zurich, Switzerland

## **Important!**

**Applying the three smaller sized B-clamps by hand may damage the calibration. Please use the clamp applying forceps for it.**

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## Clamp Boxes Klemmendosen

### Sterilization case CB-2 for original Acland®-clamps

- Easy access and clearly arranged identification of the clamps
- Suitable for all current sterilization methods
- For instant use with clamp applying forceps
- Facilitates work and classification in OR
- Provides necessary protection during sterilization, transportation and storage
- Economical, durable and functional
- Perforation provides good ventilation to facilitate drying
- A simple, cost-efficient solution to ensure longevity of your original Acland®-clamps!

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