GEM

Global Excellence in Microsurgery

NONSTICK BIPOLAR FORCEPS

Reusable Electrosurgical/Bipolar Forceps
Herbruikbare elektrochirurgische/bipolaire tang
Genanvendelig elektrokirurgisk/topolet tang
Yeniden Kullanılabilir Elektro Cerrahi/Çift Kutuplu Forsepsler
Gjenbrukbar elektrokirurgisk/bipolar tang
Forceps bipolaires électrochirurgicaux réutilisables
Wiederverwendbare elektrochirurgische/zweipolige Zange
Επαναχρησιμοποιούμενη ηλεκτροχειρουργική/διπολική λαβίδα
Pinze bipolari/per elettrochirurgia riutilizzabili
再利用型電気外科/バイポーラ鉗子
Fórceps eletrocirurgicos/bipolares reutilizáveis
Pinzas bipolares/electroquirúrgicas reutilizables
Återanvändningsbara elektrokirurgiska/bipolära peanger

CE
CE
0120

Symbol Definitions:

Attention, consult accompanying documents

Consult Instructions for Use

CONTENT

Content

Separate Collection

Rx Only

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

This product and package do not contain natural rubber latex

Non-sterile

Date of manufacture

REF

Catalog number

LOT

Lot number

EC REP

Authorized Representative in the European Community

Manufacturer

MADE IN THE U.S.A.

Made in the U.S.A.

Description

Synovis reusable bipolar forceps are designed to function as active electrodes or as accessories to active electrodes and are sold non-sterile. Electrosurgical Cords are designed to conduct electrical power from the output of a high frequency generator to a bipolar instrument. Each reusable device is designed to withstand a minimum of 20 sterilization cycles when properly cared for and sterilized according to the instructions provided below.

Indications

Bipolar Forceps and Cords are intended for use in coagulating tissue during surgical procedures and are intended for use with a maximum voltage of 500 Volts.

Handling

- Device use and application to be determined by the physician.
- Remove and discard the tip protector on all forceps prior to use.
- Inspect forceps before each use and send for repair or discard if damaged.
- Visually inspect every cord prior to each use and look for cracks, nicks, holes, or tears in the covering, as well as checking the integrity of the cord where it joins the end connectors. If any of these problems are noted or if the cord is compromised in any manner, remove the cord from service and discard.
- After the procedure, the cords should be grasped by the plug to disconnect it from the electrosurgical generator. Pulling on the cord to detach it from the generator can cause weaknesses that could lead to a fire hazard.
- Inspect bipolar forceps prior to use. Damaged electrosurgical devices may cause electrical burns to the patient or physician.
- Inspect bipolar cords for damage prior to cleaning. If nicks, holes, or tears appear in the covering of the cord or the integrity of the

connection of the cord to the plug is compromised, remove the cord from service and discard. Damaged cords could cause a fire hazard.

- Clean and decontaminate the devices according to the instructions below
 - Begin decontamination immediately after a surgical procedure.
 - Wipe off all debris with a moist sponge or gauze pad and begin decontamination immediately after the surgical procedure.
 - Use a commercially available mild detergent that is safe for plastics.
 - Cords should not be placed in an ultrasonic cleaner.
 - Rinse devices thoroughly with running tap water to remove any debris or cleaning solution.
 - Wipe dry and sterilize according to the instructions below.

Sterility

These bipolar forceps and cords are sold as non-sterile and must be steam sterilized prior to each use according to the instructions below. Other methods of sterilization are not recommended and their use will void the warranty.

Method of Steam Sterilization

Sterilization Cycle	Temperature	Time	Minimum Dry Time	Device
Prevacuum Flash	270°F (132°C)	3 Minutes	•	Sterilization Tray
Prevacuum	270°F (132°C)	4 Minutes	20 Minutes	Wrapped/ Individually Pouched
Gravity Displacement	270°F (132°C)	10 Minutes	10 Minutes	Wrapped

Storage

Reusable bipolar forceps and cords must be stored individually

in a protective container in a clean, dry, moisture free area with suitable protection from extreme temperature and humidity. Protect cords from mechanical damage and kinking and handle with care.

Cautions

- Electrosurgical instruments and accessories should be used only by individuals who have been specially trained in their use.
- Active electrosurgical devices must be placed away from patient when not in use. A holster is recommended.
- Inspect every device before each use. Damaged electrosurgical devices may cause electrical burns to the patient or physician.
- Use a wet sponge to remove eschar build up from the tips of the forceps.
- Consult RF generator Owner's Manual. Use only in the active RF generator bipolar output connections.
- Do not perform electrosurgery in the presence of flammable anesthetics or other flammable gases, fluids or objects, or in the presence of oxidizing agents, as fire could result.
- Avoid unnecessary or prolonged activation. Activate only when touching or immediately next to the intended point of contact.
- Do not bend, cut, or modify any electrosurgical device before or during surgery.
- Keep the voltage/power settings at the lowest possible level to achieve the desired effect.
- Establish a protocol for device examination and maintenance.

Warning

Improper use of any element in electrosurgical systems may cause serious burn. The use of the bipolar forceps with safely functioning and compatible cords, electrosurgical generators and accessories is the user's responsibility. Read and understand all precautions and instructions before using any electrosurgical device.

Warranty

All products are warranted to be free from defects in workmanship and materials when used for its intended purpose. Any product that proves defective in workmanship or materials will be repaired or replaced at our discretion free of charge. Repair, modification, or alteration to any product that is performed by anyone not authorized by Synovis will result in the immediate loss of warranty.

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Symmetry Surgical GmbH
Manufactured for: Synovis Micro Companies Alliance, Inc. (A subsidiary of Baxter International Inc.)