EN	NYLON	FR	NYLON
Non-Absorbable surgical suture	Fil de suture non resorbable		
Instructions For Use	Notice d'Utilisation		

DESCRIPTION: Nylon monofilament and Nylon braid sutures are inert, non-absorbable, sterile sutures composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6.6. The braided sutures are coated uniformly with silicone to enhance handling characteristics, ease of passage through tissue and reduction of capillarity. Nylon monofilament / Nylon braid sutures are non-absorbable nylon surgical sutures which are available either dyed black (with Logwood extract), undyed (clear) or dyed blue. Nylon monofilament and Nylon Braid nylon sutures meet all requirements established by the United States Pharmacopeia (USP).

INDICATIONS: Nylon monofilament and Nylon braid nylon sutures are indicated for use in general soft tissue approximation and/or ligator, including use in cardiovascular, ophthalmic, and neurological surgery. Nylon sutures are intended for one-time use only, and are not to be re-sterilized.

ACTIONS: Nylon monofilament and Nylon braid nylon sutures all elicit a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. As with any nylon suture, though not absorbent, progressive hydrolysis of the suture, in vivo, may result in gradual loss of its tensile strength over time.

CONTRAINDICATIONS: Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon sutures should not be used where permanent retention of tensile strength is required.

WARNINGS: As with any foreign body, prolonged contact of any suture whithast solutions may result in calculus formation. Do not re-sterilize. Sutures are sterile unless packaging has been opened or damaged. Discard open, unused sutures. Store at room temperature. Avoid prolonged exposure to elevated temperatures.Users should be familiar with surgical procedures and techniques involving non absorbable sutures before employing Nylon monofilament or Nylon braid sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS: In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS: Adverse effects, which may be associated with the use of this product, include: wound dehiscence, gradual loss of

DE	NYLON	IT	NYLON
Nicht resorbierbares chirurgisches Nahtmaterial	Filo di sutura chirurgico non riassorbibile		
Gebrauchsanweisung	Istruzioni per l'uso		

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tensile strength over time, calculi formation when prolonged contact with salts solutions occurs, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED: Nylon monofilament / Nylon braid nylon sutures are available either clear (undyed) or dyed black with Logwood extract, or blue, in sizes 2 (5 Metric) through 11-0 (0.1 Metric). Available in one dozen box quantities. This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or resterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize this device.

CAUTION: Federal (USA) law restricts this device to purchase and use by, or on the order of, a physician or veterinarian.

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SYMBOL DEFINITIONS	
REF	Reference Number
QTY	Quantity
STERILE EO	Sterilized by Ethylene Oxide
	Consult Instructions for Use
	Do Not Sterilize
	Do Not Reuse
	Do Not Use if Package is Damaged
	Manufacturer
LOT	Lot Number
	Expiration Date
Rx Only	

	
Riverpoint Medical 825 NE 25th Avenue Portland, OR 97232 USA Ph: +1 503-517-8001 www.rpmed.com	
	0459
	
Emergo Europe Prinsessegracht 20 2514 AP The Hague The Netherlands	

	
Global Excellence in Microsurgery	
Distributed by: Synovis Micro Companies Alliance, Inc. (A subsidiary of Baxter International Inc.) 439 Industrial Lane Birmingham, AL 35211 USA (Tel) +1 205.941.0111 (US Toll Free) 800.510.3318 (Fax) 205.941.1522 www.synovismicro.com	

PUB440-01 02-2022