NEUROCAP

Instructions for use - GB

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications

STERILE.

Sterilized with ethylene oxide gas. For single use only. Do not autoclave.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

STORAGE

- Store in a dark, dry place between -18°C and 8°C (0-46°F).
- Use the device prior to the 'use by' specified on the package.

DESCRIPTION

The NEUROCAP, a peripheral nerve capping device, is composed of the bioresorbable copolyester poly(DL-lactide-ε-caprolactone). The NEUROCAP separates the nerve from the surrounding environment and protects the nerve end.

The NEUROCAP elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the cap by fibrous tissue. Degradation of the NEUROCAP occurs through hydrolysis leading to gradual reduction of molecular weight. The NEUROCAP retains sufficient mechanical properties to act as a barrier up to at least 10 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. NEUROCAP is resorbed within approximately 16 months.

The NEUROCAP is a tubular device with one open end and one end sealed, the tip, as indicated in figure 1. An additional hole is included in the tip for easy tissue fixation with a suture. The size of the products is 3cm in length, and is available in different diameters for different nerves, see table 1.

Figure 1: NEUROCAP

The NEUROCAP inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The NEUROCAP is indicated for single-use.

Table 1: Product sizes:

Catalogue #	Length [mm]	Internal Diameter [mm]
NC01-015/03	30 ± 5	1.6 ± 0.2
NC01-020/03	30 ± 5	2.1 ± 0.2
NC01-025/03	30 ± 5	2.6 ± 0.2
NC01-030/03	30 ± 5	3.1 ± 0.3
NC01-040/03	30 ± 5	4.1 ± 0.4
NC01-050/03	30 ± 5	5.2 ± 0.4
NC01-060/03	30 ± 5	6.2 ± 0.4
NC01-070/03	30 ± 5	7.2 ± 0.4
NC01-080/03	30 ± 5	8.2 ± 0.4

INDICATIONS

The NEUROCAP is a bioresorbable device to protect a peripheral nerve end and separate the nerve from surrounding environment to reduce the development of a symptomatic neuroma.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- The NEUROCAP is for single use only. Do not re-sterilize or re-use. Structural integrity and/or function may be impaired through cleaning, re-sterilization, or re-use and may cause adverse patient

reactions. Accordingly, Polyganics will not be responsible for any direct or consequential damages or expenses resulting from re-use of (or any part of) the NEUROCAP;

- Sterile unless package has been opened or damaged. Discard open unused nerve caps;
- The NEUROCAP should only be used by surgeons who are trained in nerve defect repair techniques. Accordingly, Polyganics will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel. The physician should consult recent literature on current medical practice on peripheral nerve repair;
- Wound healing may be delayed in patients that have received radiation therapy at the site of implantation.
- Do not use the device on the volar side of the proximal interphalangeal joints

PRECAUTIONS

- Do not expose the NEUROCAP to organic solvents (e.g. chloroform, acetone);
- Do not use absorbable sutures for fixation of the nerve end into the nerve cap;
- Avoid crushing, crimping, kinking or other damage due to application of surgical instruments such as forceps, needle-holders and scissors or during handling of the device;
- Avoid tension on the nerve end;
- Ensure sufficient healthy soft tissue is available to cover the NEUROCAP in order to avoid protrusion or wound dehiscence.

ADVERSE EVENTS

Adverse events associated with the use of a NEUROCAP may include but are not limited to:

- Failure to reduce symptomatic neuroma pain
- Transitory local irritation;
- Infection;
- Allergy;
- Delayed wound healing;
- Protrusion

OPENING OF THE PACKAGE

Open the aluminum pouch. **Attention:** the outside of white/transparent pouch is not sterile. The pouch is opened in such a way that the tray remains sterile. The tray can be opened by sliding the lid. By clasping the NEUROCAP at one of its ends between a pair of tweezers, it can be taken from the tray. If it is difficult to get the device loose from the tray, soak/ rinse it in saline, preferable at 37 °C.

Inspect the cap; do not use if the device is kinked, brittle, or degraded.

NOTE: These recommendations are designed to serve only as a general procedure. They are not intended to supersede the institutional protocols or professional clinical judgment concerning patient care.

SURGICAL PROCEDURE

PRE-OPERATIVE

1. Mark the most painful spot and place of Tinel's sign before the regional block/anaesthesia has been set.

INTRA-OPERATIVE PREPARATION

1. According to hospital procedures patients receive prophylactic antibiotics.
2. Perform surgery with an upper arm tourniquet.
3. Make an incision at the marked spot.
4. Identify and dissect the nerve including neuroma.

PROCEDURE

1. Excise the neuroma with a sharp knife (the nerve stump shall NOT be diathermically transected).
2. Measure the width of the nerve.
3. Select the correct nerve cap according to the nerve diameter, see table 1. The diameter of the device should fit the diameter of the nerve. If there is no device with the size that matches the diameter of the nerve, the surgeon should select the device with 1 size up.

4. Place the NEUROCAP in warm saline (37°C) for approximately 1 minute before implantation. Temperature will make the cap more flexible and ease needle passage during suturing.
5. Cut the cap into the appropriate length with a knife or scissors up to a minimum of 15 mm, which is necessary for good positioning but leave enough length for proper positioning of the nerve (see Figure 2).
6. **For sizes 1.5 – 3.0 MM:**
7.0 or 6.0 Prolene with the smallest needle possible. Tapered Needle 3/8 (9 – 11 mm).
For sizes 4.0-8.0 MM:
5.0 or 6.0 – Prolene or monofilament with 11 mm tapered needle.
5.0 or 6.0 Ethilon with 13 mm needle or with the smallest tapered needle available.
7. Start the first stitch from the outside of the nerve cap, 10 mm from the cap tip extremity.
8. Take the epineurium at 5mm of the nerve end and then back from inside the nerve cap, outside.
9. If it is not possible in the surgeon's opinion to glide the nerve into the cap easily with one suture, repeat the same steps on the other side of the cap and nerve.
10. After this position the nerve is surrounded with sufficient soft tissue, preferably muscle.
11. Use the tip end (with hole) for fixating the nerve cap ending by using a 5-0 ethilon or prolene suture onto the soft tissue, fascia, muscle or periost.


Figure 2. Schematic representation of suture technique for suturing the nerve end into the nerve cap.

POST-OPERATIVE

The patients should receive a bandage dressing. If movement of the nerve, due to joint movement, causes tension on the device, a plaster might be indicated. If the NEUROCAP was placed over a joint a plaster might be indicated. Advice to patients not to put any pressure on the surgery side, to avoid compression or dislocation of the nerve during 6 weeks postoperative. In case of splinting, splints also should be on the other side of the scar. The patient receives a sling and adequate painkillers for postoperative pain relief.

DISPOSAL

Dispose contaminated implantation and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

The NEUROCAP is MR Safe.

USED SYMBOLS

- = Consult instructions for use.
- = Do not use if package is damaged
- = Manufacturer
- = Use by
- REF** = Catalog number
- = Do not reuse
- = Do not re-sterilize
- LOT** = Lot number
- STERILE EO** = Sterilized with ethylene oxide gas
- = Keep dry
- = Keep away from sunlight
- = Storage temperature limit
- MR** = MR Safe

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

- There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Polyganics product(s) described in this publication.
- Under no circumstances shall Polyganics be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law.
- No person has the authority to bind Polyganics to any representation or warranty except as specifically set forth herein.
- Descriptions or specifications in Polyganics printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.
- Polyganics is not responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.

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